


BFI USB style
Cell Check Plus II



【PC windows10】




BFI Software USB
Probe USB style 7.5MHz



Model number	SM706
Product name	BFI Measure Cellu Check Plus (Ultrasound imagemeasremnt device)
Construction	USB Probe/BFI USB software Attachment parts, Sonic water/User manual, Certification, Check sheet
Measurement parts	Upper arm, Abdomen, Waist, Front thigh, Back thigh
Measurement information	1 millimeter (mm), Water retention (%)
Communication method	USB
Measurements items	Fat thickness, Muscle thickness, Cell check, FAT total, FAT%, FAT kg, Body Mass Index, Muscle weight, Obesity index, Basal metabolism, Body Balance, Metabolic, Diet plan, Body balance, Delete calorie plan, Food & Exercise plan, Target BML, Target simulation, Image simulation, Save measurement data, Progress data, Check sheet (Print sheet)
Result check sheet	Special paper (check sheet A4size)
Measurement sound	Touch panel sound attached * It could be changed by PC operation.
Result Save	Image, a numerical saved
Measurement method	Pulse method Result Save Image, a numerical saved
Probe	7.5Mhz 128elements linear model
Measurement wide	About 75mm
Measurement depth	About 70mm
Size	Probe 102mm(W) × 40mm(D) × 47(H)mm
Weight	Probe 300g(include cable)
Soft ware	Body Fat Index software
Product name	BFI706
Method	USB memory * The management software is included in normal constitution.

Soft ware Language 【Japan, English, French, Spanish】



FAMUBO



Size	60(W) × 24(H) × 156(L) mm Weight: 0.3kg
Battery	3 hour
Probe	7.5MHz / 10MHz 128 elements 32channel
Probe width size	Linear 40mm
Depth	linear 40- 100mm / Display Mod B
OS spec	iPad mini 4 ~ / iPad air 2 ~ / iPad Pro ~
OS spec	iPhone 5 ~ / Android 6.0 ~

Soft ware Language 【Japan, English】



FAMUBO



BFI USB style
Cell Check Plus II

SEIKOSHA
Innovative Republic

New Body composition meter by Ultrasound

Seikosha Co., LTD. studied ultrasound from 1990 and analyzed composition of 8cm under the skin and announced 2 kinds the new body composition analysis by the latest technology and abundant software.

"Anytime" "Anywhere" "Anyone"
It's very easy!!
It's reasonable price!!
It's very high quality picture image!
It's a compact!!



PTDR. NAGAKI




BFI USB style
Cell Check Plus II

BFI (Body Fat Index) New index for Body

The BMI is calculate from weight and height.
(BMI=Weight (kg) / height (m)²)
But if there is same BMI people, one people has big abdomen, the other people has big arms and legs. It is related their lifestyle. BMI is not reflected body fat at all. Therefore, it could not possible to advice the counseling with body balance details. The weight and Fat percentage are insufficiency for details counseling. It find just obesity level and it couldn't find body balance for part of body. The BFI is added fat thickness with weight and height by ultrasound. Where the fat thickness is most big? How is the body balance? The BFI could find the figure and body balance and could advice details date of customer.

- ①BFI Balance ②Fat % ③FAT Thickness & Orange Peel ④Diet Calorie Plan will be able to measured as a measurement items. It is possible to management customer, patient, client and analyzes and evaluates the 6 blocks balance of Body Figure balance, Body's figure evaluation, Part of Fat Thickness balance, Diet Calorie Plan.



Anti-aging Dr. AOII

FAMUBO (Wireless Ultrasound)

Ultrasound image measurement has evolved !!
The "FAMUBO", has appeared as small wireless type(Wi-Fi) and reasonable price. "Every time", "anywhere", "anyone" can use quickly and easily. The image quality is very high and it's as good as the image quality of large and expensive models. Amazingly, "FAMUBO" can use as the body composition meter too. By simply connecting to your existing IPAD, it can be quickly turned into an ultrasound image measurement device and Body composition meter.

You are Wireless style or USB style?

Measurement Item

FAMUBO have a function both ultrasound device and body composition meter.
①Ultrasound image ②Fat % FAT (body composition meter)
③Diet Calorie Plan

Registration

Name, SEX, address, TEL, email, Age, height(cm)
weight (kg), waist (cm)

Ultrasound Image

The picture is very high quality and it's a almost same level compare with expensive device.

Depth: 4cm, 6cm, 8cm, 10cm
Probe: 7.5MHz & 10MHz
Image : Black & white image




Measurement Item

①BFI Balance ②Fat % ③FAT Thickness & Orange Peel ④Diet Calorie Plan It is possible to measure four items measurement.

Registration

Age, height (cm), weight (kg), waist (cm)



Body's Figure Balance

After calculation, it will be indicated the Result of 6 blocks of FAT% and BMI. A vertical axis is FAT%, a cross axis is BMI. You could find the balance on this screen.

Part of Fat Thickness balance

You could find the part of Fat thickness balance (BFI-part of BMI) on this list. It indicates that where is your fattest part and where is your thinnest part and evaluates BFI (part of BMI). This table is best for diet program and therapy method for body part.

Result sheet

FAT %, FAT kg, Muscle kg, BMI, Basal metabolic rate
FAT thickness (free) Muscle thickness (free)
Doctor comment

DATE save

You can compare with 2 picture image on IPAD.
You can overlap two ultrasound image.
You can compare on left and right.

Progress Data (Check my data)

It is saved your date of every measurement. If you click the measurement data, you could find the ultrasound images which you measured

Simulation image

If you click the "SIMULATIO" button, You could find the future slim body image in a moment. If you change the BMI, you could find the more slim image or more big image. You can keep the motivation of customer..

Diet Calorie Plan

You could calculate calorie plan. At first, please input "TARGET WEIGHT", after that, please input term, you could find the calorie per day. If you input the calorie, you could find the term. You can enclose your customer for a few months.

Ultrasound Color analysis Image

It project a ultrasound image on PC and could see subcutaneous fat, muscle, cell fat, bone with ultrasound visual. You can judge sports type or cellulite type or soft type or obesity type for each parts for body.

Body's Figure Balance

After calculation, it will be indicated the Result of 6 blocks of FAT% and BMI. A vertical axis is FAT%, a cross axis is BMI. You could find the balance on this screen. The accuracy is more than 90%.

