# ESTHELUX-Pro

PHOTO BIO MODULATION by Doctor

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

Dr. In Seung-gyun  
Human Dermatology

**Effect of Light Emitting Diode Photomodulation in Reducing Erythema After Fractional Carbon Dioxide Laser Resurfacing: A Pilot Study**

Background: Fractional carbon dioxide laser resurfacing is a common procedure for skin rejuvenation. However, it often causes post-inflammatory erythema (PIE), which can be a significant concern for patients. Light-emitting diode (LED) photomodulation has been shown to have anti-inflammatory and vascular effects. This study aimed to evaluate the effect of LED photomodulation on reducing erythema after fractional carbon dioxide laser resurfacing.

Methods: A pilot study was conducted with 10 patients who underwent fractional carbon dioxide laser resurfacing. The patients were randomized into two groups: one group received LED photomodulation immediately after the laser treatment, and the other group did not receive LED photomodulation. The erythema was assessed using a visual analog scale (VAS) at baseline, immediately after the procedure, and at 1 hour, 24 hours, and 48 hours post-procedure.

Results: The LED photomodulation group showed significantly lower VAS scores compared to the control group at all time points. The difference was most pronounced at 1 hour and 24 hours post-procedure.

Conclusion: LED photomodulation effectively reduces erythema after fractional carbon dioxide laser resurfacing. This finding suggests that LED photomodulation may be a useful adjunct to laser resurfacing to improve patient outcomes and reduce recovery time.

Dr. Jeong Chan-woo  
JF Dermatology

**Correspondence**


After photomodulation with LED, the erythema was significantly reduced. This finding suggests that LED photomodulation may be a useful adjunct to laser resurfacing to improve patient outcomes and reduce recovery time.

Dr. Kim Ki-bum  
Misogain Dermatology

## A case of small plaque psoriasis Treated with Combination of 578-nm Copper Bromide Laser with Light-emitting Diode (LED)

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